



MEDSTIM

QuickFit TTFM Probes

Verified quality.
Better outcomes.



Reliable, proven technology for improving surgical outcome

Medistim's flow probes are based on reliable Transit Time Flow Measurement (TTFM). Together with the VeriQ™ and VeriQ C™ systems, they provide fast, accurate and reproducible measurements of blood volume flow intraoperatively. This is particularly relevant to verifying graft patency and function in a range of surgical applications: cardiac, vascular and transplants. The ultimate benefit is quality assurance with improved surgical outcome.

Specifically designed for cardiac surgery

You are able to measure all vessels 1.3–6.0 mm in size. Probes are available with and without handles in order to accommodate a variety of surgical procedures. Designed to meet worldwide sterilization standards, Medistim's QuickFit probes are sterilizable by gas (EtO) or plasma (STERRAD®) and are guaranteed for 50 sterilization cycles. The probes are factory calibrated and require no further calibration during use.

QuickFit probes with handle

Probe size mm	Size of vessel/graft mm	System Compatibility	Part number
1.5	1.3–2.0	VeriQ C, VeriQ, ButterFly	PQ101012
2	1.8–3.0	VeriQ C, VeriQ, ButterFly	PQ100022
3	2.4–4.0	VeriQ C, VeriQ, ButterFly	PQ100032
4	3.2–5.0	VeriQ C, VeriQ, ButterFly	PQ100042
5	4.0–6.0	VeriQ C, VeriQ, ButterFly	PQ100052

QuickFit probes without handle

Probe size mm	Size of vessel/graft mm	System Compatibility	Part number
1.5	1.3–2.0	VeriQ C, VeriQ, ButterFly	PQ101011
2	1.8–3.0	VeriQ C, VeriQ, ButterFly	PQ100021
3	2.4–4.0	VeriQ C, VeriQ, ButterFly	PQ100031
4	3.2–5.0	VeriQ C, VeriQ, ButterFly	PQ100041
5	4.0–6.0	VeriQ C, VeriQ, ButterFly	PQ100051

Please refer to the User Manual for indications, contraindications, warnings, precautions, and further specifications and descriptions. Specifications may be changed without notice. For a list of flow probes for other applications, contact your Medistim representative.

STERRAD®

FDA 510(k) cleared no. K102595



medistim@medistim.com
www.medistim.com

Medistim ASA (Head office)
Økernveien 94
0579 Oslo
Norway
Phone +47 23 05 96 60

Medistim Norge AS
Økernveien 94
0579 Oslo
Norway
Phone +47 23 03 52 50

Medistim USA Inc.
14000 25th Ave N. Ste. 108
Plymouth, MN 55447
USA
Phone +1 763 208 9852

Medistim Deutschland GmbH
Bahnhofstr. 32
82041 Deisenhofen
Germany
Phone +49 (0) 89 62 81 90 33

Medistim Danmark ApS
Gøngetoften 13
2950 Vedbæk
Denmark
Phone +45 2276 5669

Medistim UK Limited
2nd Floor
2 City Place Beehive Ring Road
Gatwick
West Sussex RH6 0PA
Phone: +44 207 486 5888