

# Vascular TTFM Probes

Verified quality.  
Better outcomes.

**Reliable, proven technology for improving surgical outcome**

Medistim's flow probes are based on reliable Transit Time Flow Measurement (TTFM). Together with the VeriQ™ and VeriQ C™ systems, they provide fast, accurate and reproducible measurements of blood volume flow intraoperatively. This is particularly relevant to verifying graft patency and function in a range of surgical applications: cardiac, vascular and transplants. The ultimate benefit is quality assurance with improved surgical outcome.

**Specifically designed for peripheral vascular and transplant surgery**

You are able to measure all vessels 5–18 mm in size. Probes are available with and without handles in order to accommodate a variety of surgical procedures. Designed to meet worldwide sterilization standards, Medistim's vascular probes are sterilizable by gas (EtO) or plasma (STERRAD®) and are guaranteed for 50 sterilization cycles. The probes are factory calibrated and require no further calibration during use.

**Vascular probes with handle**

Probe size mm	Size of vessel/graft mm	System Compatibility	Part number
6	5–7	VeriQ C, VeriQ, ButterFly	PB100062
8	6–9	VeriQ C, VeriQ, ButterFly	PA100082
10	8–12	VeriQ C, VeriQ, ButterFly	PA100102
12	9–13	VeriQ C, VeriQ, ButterFly	PA100122

**Vascular probes without handle**

Probe size mm	Size of vessel/graft mm	System Compatibility	Part number
8	6–9	VeriQ C, VeriQ, ButterFly	PA100081
10	8–12	VeriQ C, VeriQ, ButterFly	PA100101
12	9–13	VeriQ C, VeriQ, ButterFly	PA100121
16	12–18	VeriQ C, VeriQ, ButterFly	PA100161

Please refer to the User Manual for indications, contraindications, warnings, precautions, and further specifications and descriptions. Specifications may be changed without notice. For a list of flow probes for other applications, contact your Medistim representative.

STERRAD®

FDA 510(k) cleared no. K102595



[medistim@medistim.com](mailto:medistim@medistim.com)  
[www.medistim.com](http://www.medistim.com)

**Medistim ASA (Head office)**  
 Økernveien 94  
 0579 Oslo  
 Norway  
 Phone +47 23 05 96 60

**Medistim Norge AS**  
 Økernveien 94  
 0579 Oslo  
 Norway  
 Phone +47 23 03 52 50

**Medistim Danmark ApS**  
 Gøngetoften 13  
 2950 Vedbæk  
 Denmark  
 Phone +45 2276 5669

**Medistim USA Inc.**  
 14000 25th Ave N. Ste. 108  
 Plymouth, MN 55447  
 USA  
 Phone +1 763 208 9852

**Medistim Deutschland GmbH**  
 Bahnhofstr. 32  
 82041 Deisenhofen  
 Germany  
 Phone +49 (0) 89 62 81 90 33

**Medistim UK Limited**  
 2 City Place  
 Beehive Ring Road, Gatwick  
 West Sussex, England  
 Phone +44 207 486 5888